

Product datasheet for **TP306392M**

FGFBP2 (NM_031950) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human fibroblast growth factor binding protein 2 (FGFBP2), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206392 protein sequence Red =Cloning site Green =Tags(s) MKFVPCLLLVTLSCLGTLGQAPRQKQGSTGEEFHFQTGGRDCTMRPSSLGQGAGEVWLRVDCRNTDQTY WCEYRGQPSMCQAFADPKPYWNQALQELRRLHHACQGAPVLRPSVCREAGPQAHMQVTSLLKGSPEPN QQPEAGTPSLRPKATVKLLEATQLGKDSMEELGKAKPTTRPTAKPTQPGPRPGNNEAKKKAWHCWKPF QALCAFLISFFRG TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	22.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_114156
Locus ID:	83888
UniProt ID:	Q9BYJ0



[View online »](#)

RefSeq Size: 1188

Cytogenetics: 4p15.32


RefSeq ORF: 669

Synonyms: HBP17RP; KSP37

Summary: This gene encodes a member of the fibroblast growth factor binding protein family. The encoded protein is a serum protein that is selectively secreted by cytotoxic lymphocytes and may be involved in cytotoxic lymphocyte-mediated immunity. An increase in the amount of gene product may be associated with atopic asthma and mild extrinsic asthma.[provided by RefSeq Staff, Oct 2008]

Protein Families: Secreted Protein

Product images:



Coomassie blue staining of purified FGFBP2 protein (Cat# [TP306392]). The protein was produced from HEK293T cells transfected with FGFBP2 cDNA clone (Cat# [RC206392]) using MegaTran 2.0 (Cat# [TT210002]).