

## Product datasheet for **TP306402**

### HDHD3 (NM\_031219) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human haloacid dehalogenase-like hydrolase domain containing 3 (HDHD3), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206402 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	 MAHRLQIRLLTWDVKDTLLRLRHPLGEAYATKARAHGLEVEPSALEQGFRQAYRAQSHSFPNYGLSHGLT SRQWWLDVWLQTFHLAGVQDAQAVAPIAEQLYKDFSHPCTWQVLDGAEDTLRECRTRGLRLAVISNFDRR LEGILGGLGLREHDFVLTSEAAGWPKDPDRIFQEALRLAHMEPVVAAHVGDNYLCDYQGPRAVGMHSFL VWGPQALDPVVRDSVPKEHILPSLAHLLPALDCLEGSTPGL  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	27.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_112496</a></u>
Locus ID:	81932



[View online »](#)

UniProt ID: [Q9BSH5](#)  
RefSeq Size: 1787  
Cytogenetics: 9q32  
RefSeq ORF: 753  
Synonyms: 2810435D12Rik; C9orf158

**Product images:**

Coomassie blue staining of purified HDHD3 protein (Cat# TP306402). The protein was produced from HEK293T cells transfected with HDHD3 cDNA clone (Cat# [RC206402]) using MegaTran 2.0 (Cat# [TT210002]).