

Product datasheet for **TP306411L**

GRRP1 (FAM110D) (NM_024869) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human glycine/arginine rich protein 1 (GRRP1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206411 protein sequence Red =Cloning site Green =Tags(s)

MLLAPPSTPSRGRTPSAVERLEADKAKYVKTHQVIARRQEPALRGSPGLTPHPCNELGPPASPRTPRPV
RRGSGRRLPRPDSLIFYRQKRDCKASVKNENAKGQGLVRRFLGAPRDAAPSSPASTERPAASGGWAAPQ
DAPEAAGKRALCPTCSLPLSEKERFFNYCGLERALVEVLGAERFSPQSWGADASPQAGTSPPPGSGDASD
WTSSDRGVDSPPGGAGGGGGSEAAGSARDRRPPVSVVERNARVIQWLYGCQRARGPPRESEV


TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	28.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_079145
Locus ID:	79927
UniProt ID:	Q8TAY7



[View online »](#)

RefSeq Size: 1568
Cytogenetics: 1p36.11
RefSeq ORF: 813
Synonyms: GRRP1

Product images:

Coomassie blue staining of purified FAM110D protein (Cat# [TP306411]). The protein was produced from HEK293T cells transfected with FAM110D cDNA clone (Cat# [RC206411]) using MegaTran 2.0 (Cat# [TT210002]).