

## Product datasheet for **TP306479L**

### TRIM (TRAT1) (NM\_016388) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human T cell receptor associated transmembrane adaptor 1 (TRAT1), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206479 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)  MSGISGCPFFLWGLLALLGLALVISLIFNISHYVEKQRQDKMYSYSSDHTRVDEYYIEDTPIYGNLDDMI SEPMDENCYEQMKARPEKSVNKMQEATPSAQATNETQMCYASLDHSVKGKRRKPRKQNTHFSDKDGDEQL HAIDASVSKTTLVDSFSPESQAVEENIHDDPIRLFGLIRAKREPIN  <b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	21 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_057472</a>
Locus ID:	50852
UniProt ID:	<a href="#">Q6PIZ9</a>



[View online »](#)

RefSeq Size:	1695
Cytogenetics:	3q13.13
RefSeq ORF:	558
Synonyms:	HSPC062; pp29/30; TCRIM; TRIM
Summary:	Stabilizes the TCR (T-cell antigen receptor)/CD3 complex at the surface of T-cells. [UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome, Transmembrane

**Product images:**

Coomassie blue staining of purified TRAT1 protein (Cat# [TP306479]). The protein was produced from HEK293T cells transfected with TRAT1 cDNA clone (Cat# [RC206479]) using MegaTran 2.0 (Cat# [TT210002]).