

Product datasheet for **TP306558L**

HSD17B3 (NM_000197) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hydroxysteroid (17-beta) dehydrogenase 3 (HSD17B3), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206558 protein sequence Red =Cloning site Green =Tags(s)

MGDVLEQFFILTGLLVCLACLAKCVRFRCVLLNYWKVLPKSFRLSMGQWAVITGAGDGIGKAYSFELAK
RGLNVWLISRTLEKLEAIAIEIERTTGRSVKIIQADFTKDDIYEHKEKLAGLEIGILVNNVGMPLNLLP
SHFLNAPDEIQSLIHCNITSVVKMTQLILKHMESRQKGLILNISSGIALFPWPLYSMYSASKAFVCAFSK
ALQEEYKAKEVIIQVLTPYAVSTAMTKYLNTNVITKTADEFVKESLNYVTIGGETCGCLAHEILAGFLSL
IPAWAFYSGAFQRLLLTHYVAYLKLNTKVR

TRTRPLE**QKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	34.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_000188</u>
Locus ID:	3293



[View online »](#)

UniProt ID: [P37058](#), [Q6FH62](#), [G5E9S2](#)

RefSeq Size: 1134

Cytogenetics: 9q22.32

RefSeq ORF: 930


Synonyms: EDH17B3; SDR12C2

Summary: This isoform of 17 beta-hydroxysteroid dehydrogenase is expressed predominantly in the testis and catalyzes the conversion of androstenedione to testosterone. It preferentially uses NADP as cofactor. Deficiency can result in male pseudohermaphroditism with gynecomastia. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transmembrane

Protein Pathways: Androgen and estrogen metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified HSD17B3 protein (Cat# [TP306558]). The protein was produced from HEK293T cells transfected with HSD17B3 cDNA clone (Cat# [RC206558]) using MegaTran 2.0 (Cat# [TT210002]).