

## Product datasheet for **TP306562**

### CD247 (NM\_000734) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human CD247 molecule (CD247), transcript variant 2, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206562 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MKWKALFTAAILQAQLPITEAQSFGLLDPKLCYLLDGILFIYGVILTALFLRVKFSRSADAPAYQQGQNNQ LYNELNLGRREEYDVLDKRRGRDPEMGGKPKRRKPNQEGLYNELQKDKMAEAYSEIGMKGERRRGKGGH DG LYQGLSTATKDTYDALHMQALPPR
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	16.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u><a href="#">NP_000725</a></u>
Locus ID:	919
UniProt ID:	<u><a href="#">P20963</a></u>



[View online »](#)

RefSeq Size: 1687

Cytogenetics: 1q24.2

RefSeq ORF: 492


Synonyms: CD3-ZETA; CD3H; CD3Q; CD3Z; IMD25; T3Z; TCRZ

**Summary:** The protein encoded by this gene is T-cell receptor zeta, which together with T-cell receptor alpha/beta and gamma/delta heterodimers, and with CD3-gamma, -delta and -epsilon, forms the T-cell receptor-CD3 complex. The zeta chain plays an important role in coupling antigen recognition to several intracellular signal-transduction pathways. Low expression of the antigen results in impaired immune response. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome, Transmembrane

**Protein Pathways:** Natural killer cell mediated cytotoxicity, T cell receptor signaling pathway

### Product images:



Coomassie blue staining of purified CD247 protein (Cat# TP306562). The protein was produced from HEK293T cells transfected with CD247 cDNA clone (Cat# [RC206562]) using MegaTran 2.0 (Cat# [TT210002]).