

## Product datasheet for **TP306564M**

### ASPA (NM\_000049) Human Recombinant Protein

#### Product data:

|                                       |  |
|---------------------------------------|--|
| Product Type:                         | Recombinant Proteins   |
| Description:                          | Recombinant protein of human aspartoacylase (Canavan disease) (ASPA), transcript variant 1, 100 µg |
| Species:                              | Human  |
| Expression Host:                      | HEK293T  |
| Expression cDNA Clone or AA Sequence: | >RC206564 protein sequence<br><b>Red</b> =Cloning site <b>Green</b> =Tags(s)                       |

MTSCHIAEEHIQKVAIFGGTHGNELTGVFLVKHWLENGAEIQRGTGLEVKPFITNPRAVKKCTRYIDCDLN  
RIFDLENLGKKMSEDLPEYVRRRAQEINHLFGPKDSEDSYDIIFDLHNTTSNMGCTLILEDNRNNFLIQMF  
HYIKTSLAPLPCYVYLIEHPSLKYATTRSIAKYPVGIEVGPQPQGVLRADILDQMRMKIKHALDFIHHFN  
EGKEFPPCAIEVYKIIKVDYPRDENGIEAIIHPNLQDQDWKPLHPGDPMFLTLDGKTIPLGGDCTVYP  
VFNAAAYYEKKEAFKTTKLTNNAKSIRCCLH

**TR**TRPLEQKLISEEDLAANDILDYKDDDDKV

|                |  |
|----------------|--|
| Tag:           | C-Myc/DDK  |
| Predicted MW:  | 35.6 kDa   |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method   |
| Purity:        | > 80% as determined by SDS-PAGE and Coomassie blue staining  |
| Buffer:        | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol   |
| Preparation:   | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.                                     |
| Note:          | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage:       | Store at -80°C.  |
| Stability:     | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.        |
| RefSeq:        | <u><a href="#">NP_000040</a></u>   |
| Locus ID:      | 443  |



[View online »](#)

UniProt ID: [P45381, Q6FH48](#)

RefSeq Size: 1435

Cytogenetics: 17p13.2

RefSeq ORF: 939


Synonyms: ACY2; ASP

**Summary:** This gene encodes an enzyme that catalyzes the conversion of N-acetyl\_L-aspartic acid (NAA) to aspartate and acetate. NAA is abundant in the brain where hydrolysis by aspartoacylase is thought to help maintain white matter. This protein is an NAA scavenger in other tissues. Mutations in this gene cause Canavan disease. Alternatively spliced transcript variants have been found for this gene. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome

**Protein Pathways:** Alanine, aspartate and glutamate metabolism, Histidine metabolism

### Product images:



Coomassie blue staining of purified ASPA protein (Cat# [TP306564]). The protein was produced from HEK293T cells transfected with ASPA cDNA clone (Cat# [RC206564]) using MegaTran 2.0 (Cat# [TT210002]).