

Product datasheet for **TP306565L**

ARSB (NM_198709) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human arylsulfatase B (ARSB), transcript variant 2, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206565 protein sequence Red =Cloning site Green =Tags(s)
	<p>MGPRGAASLPRGPGPRLLLLPVVLPLLLLLLAPPGSGAGASRPPHLVFLADDLGWNDVGFHGSRIPT HLDALAAGGVLLDNYYTQPLCTPSRSQLLTGRYQIRTGLQHQQIWPQPSCVPLDEKLLPQLLKEAGYTT HMVGKWHLGMRYRKECLPTRRGFDYFGYLLGSEDYYSHERCTLIDALNVTRCALDFRDGEEVATGYKNMY STNIFTKRAIALITNHPPEKPLFLYLALQSVHEPLQVPEEYLKPYDFIQDKNRHHYAGMVSLMDEAVGNV TAALKSSGLWNNTVFIFSTDNGGQTLAGNNWPLRGRKWSLWEGGVRGVGFVAPLLKQKGVKNRELIHI SDWLPTLVKLGARGHTNGTKPLDGFVWKTISEGSPSPRIELLHNIDPNFVDSSPYWPECSLLL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	42.3 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_942002</u>
Locus ID:	411



[View online »](#)

UniProt ID: [P15848](#), [A8K4A0](#)

RefSeq Size: 2172

Cytogenetics: 5q14.1

RefSeq ORF: 1239


Synonyms: ASB; G4S; MPS6

Summary: Arylsulfatase B encoded by this gene belongs to the sulfatase family. The arylsulfatase B homodimer hydrolyzes sulfate groups of N-Acetyl-D-galactosamine, chondroitin sulfate, and dermatan sulfate. The protein is targeted to the lysosome. Mucopolysaccharidosis type VI is an autosomal recessive lysosomal storage disorder resulting from a deficiency of arylsulfatase B. Two alternatively spliced transcript variants encoding distinct isoforms have been found for this gene. [provided by RefSeq, Dec 2016]

Protein Families: Druggable Genome

Protein Pathways: Glycosaminoglycan degradation, Lysosome, Metabolic pathways

Product images:



Coomassie blue staining of purified ARSB protein (Cat# [TP306565]). The protein was produced from HEK293T cells transfected with ARSB cDNA clone (Cat# [RC206565]) using MegaTran 2.0 (Cat# [TT210002]).