

## Product datasheet for **TP306572M**

### Aurora A (AURKA) (NM\_198436) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human aurora kinase A (AURKA), transcript variant 5, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206572 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)

MDRSKENCISGPVKATAPVGGPKRVLVTQQFPCQNPLPVNSGQAQRVLCPSNSSQRVPLQAQKLVSSHKP  
VQNQKQKQLQATSVPHVSRPLNNTQKSKQPLPSAPENNP EEELASKQKNEESKKRQWALEDFEIGRPLG  
KKGKFGNVYLAREKQSKFILALKVLFKAQLEKAGVEHQLRREVEIQSHLRHPNILRLYGYFHDATRVYLIL  
EYAPLGTVYRELQKLSKFDEQRTATYITELANALSYCHSKRVIHRDIKPENLLLGSAGELKIADFGWSVH  
APSSRRTTLCGTLDYLPPEMIEGRMHDEKVDLWLSGLVLCYEFLVGKPPFEANTYQETYKRISRVFTFPD  
FVTEGARDLISRLKHNPSQRPMLREVLEHPWITANSSKPSNCQNKESASKQS

**TRTRPLEQKLISEEDLAANDILDYKDDDDKV**

Tag:	C-Myc/DDK
Predicted MW:	45.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_940838</a>
Locus ID:	6790



[View online »](#)

UniProt ID: [Q14965](#)

RefSeq Size: 2231

Cytogenetics: 20q13.2

RefSeq ORF: 1209


Synonyms: AIK; ARK1; AURA; BTAK; PPP1R47; STK6; STK7; STK15

**Summary:** The protein encoded by this gene is a cell cycle-regulated kinase that appears to be involved in microtubule formation and/or stabilization at the spindle pole during chromosome segregation. The encoded protein is found at the centrosome in interphase cells and at the spindle poles in mitosis. This gene may play a role in tumor development and progression. A processed pseudogene of this gene has been found on chromosome 1, and an unprocessed pseudogene has been found on chromosome 10. Multiple transcript variants encoding the same protein have been found for this gene. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome, Protein Kinase, Stem cell - Pluripotency

**Protein Pathways:** Oocyte meiosis

### Product images:



Coomassie blue staining of purified AURKA protein (Cat# [TP306572]). The protein was produced from HEK293T cells transfected with AURKA cDNA clone (Cat# [RC206572]) using MegaTran 2.0 (Cat# [TT210002]).