

Product datasheet for **TP306581**

RBMS2 (NM_002898) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human RNA binding motif, single stranded interacting protein 2 (RBMS2), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206581 protein sequence Red =Cloning site Green =Tags(s) <p> M L L S V T S R P G I S T F G Y N R N N K K P Y V S L A Q Q M A P P S P S N S T P N S S S G S N G N D Q L S K T N L Y I R G L Q P G T T D Q D L V K L C Q P Y G K I V S T K A I L D K T T N K C K G Y G F V D F D S P S A A Q K A V T A L K A S G V Q A Q M A K Q Q E Q D P T N L Y I S N L P L S M D E Q E L E G M L K P F G Q V I S T R I L R D T S G T S R G V G F A R M E S T E K C E A I I T H F N G K Y I K T P P G V P A P S D P L L C K F A D G G P K K R Q N Q G K F V Q N G R A W P R N A D M G V M A L T Y D P T T A L Q N G F Y P A P Y N I T P N R M L A Q S A L S P Y L S S P V S S Y Q R V T Q T S P L Q V P N P S W M H H S Y L M Q P S G S V L T P G M D H P I S L Q P A S M M G P L T Q Q L G H L S L S S T G T Y M P T A A A M Q G A Y I S Q Y T P V P S S S V S V E E S S G Q Q N Q V A V D A P S E H G V Y S F Q F N K T R T R P L E Q K L I S E E D L A A N D I L D Y K D D D D K V </p>
Tag:	C-Myc/DDK
Predicted MW:	43.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_002889</u>
Locus ID:	5939
UniProt ID:	<u>Q15434</u>
RefSeq Size:	8504
Cytogenetics:	12q13.3
RefSeq ORF:	1221
Synonyms:	SCR3
Summary:	The protein encoded by this gene is a member of a small family of proteins which bind single stranded DNA/RNA. These proteins are characterized by the presence of two sets of ribonucleoprotein consensus sequence (RNP-CS) that contain conserved motifs, RNP1 and RNP2, originally described in RNA binding proteins, and required for DNA binding. The RBMS proteins have been implicated in such diverse functions as DNA replication, gene transcription, cell cycle progression and apoptosis. This protein was isolated by phenotypic complementation of cdc2 and cdc13 mutants of yeast and is thought to suppress cdc2 and cdc13 mutants through the induction of translation of cdc2. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified RBMS2 protein (Cat# TP306581). The protein was produced from HEK293T cells transfected with RBMS2 cDNA clone (Cat# [RC206581]) using MegaTran 2.0 (Cat# [TT210002]).