

Product datasheet for TP306621L

OriGene Technologies, Inc.

9620 Medical Center Drive, Ste 200 Rockville, MD 20850, US Phone: +1-888-267-4436 https://www.origene.com techsupport@origene.com EU: info-de@origene.com CN: techsupport@origene.cn

ALOX15 (NM 001140) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human arachidonate 15-lipoxygenase (ALOX15), 1 mg

Species: Human
Expression Host: HEK293T

Expression cDNA Clone >RC206621 representing NM_001140 or AA Sequence: Red=Cloning site Green=Tags(s)

MGLYRIRVSTGASLYAGSNNQVQLWLVGQHGEAALGKRLWPARGKETELKVEVPEYLGPLLFVKLRKRHL LKDDAWFCNWISVQGPGAGDEVRFPCYRWVEGNGVLSLPEGTGRTVGEDPQGLFQKHREEELEERRKLYR WGNWKDGLILNMAGAKLYDLPVDERFLEDKRVDFEVSLAKGLADLAIKDSLNVLTCWKDLDDFNRIFWCG QSKLAERVRDSWKEDALFGYQFLNGANPVVLRRSAHLPARLVFPPGMEELQAQLEKELEGGTLFEADFSL LDGIKANVILCSQQHLAAPLVMLKLQPDGKLLPMVIQLQLPRTGSPPPPLFLPTDPPMAWLLAKCWVRSS DFQLHELQSHLLRGHLMAEVIVVATMRCLPSIHPIFKLIIPHLRYTLEINVRARTGLVSDMGIFDQIMST GGGGHVQLLKQAGAFLTYSSFCPPDDLADRGLLGVKSSFYPQDALRLWEIIYRYVEGIVSLHYKTDVAVK DDPELQTWCREITEIGLQGAQDRGFPVSLQARDQVCHFVTMCIFTCTGQHASVHLGQLDWYSWVPNAPCT MRLPPPTTKDATLETVMATLPNFHQASLQMSITWQLGRRQPVMVAVGQHEEEYFSGPEPKAVLKKFREEL

AALDKEIEIRNAKLDMPYEYLRPSVVENSVAI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 74.6 kDa

Concentration: >0.05 μg/μL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional

chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience

some loss of protein during the filtration process.

Storage: Store at -80°C.





ALOX15 (NM_001140) Human Recombinant Protein - TP306621L

Stability: Stable for 12 months from the date of receipt of the product under proper storage and

handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: NP 001131

Locus ID: 246

 UniProt ID:
 P16050

 RefSeq Size:
 2707

Cytogenetics: 17p13.2 RefSeq ORF: 1986

Synonyms: 12-LOX; 15-LOX; 15-LOX-1; LOG15

Summary: This gene encodes a member of the lipoxygenase family of proteins. The encoded enzyme acts


on various polyunsaturated fatty acid substrates to generate various bioactive lipid mediators such as eicosanoids, hepoxilins, lipoxins, and other molecules. The encoded enzyme and its reaction products have been shown to regulate inflammation and immunity. Multiple pseudogenes of this gene have been identified in the human genome. [provided by RefSeq,

Aug 2017]

Protein Families: Druggable Genome

Protein Pathways: Arachidonic acid metabolism, Linoleic acid metabolism, Metabolic pathways

Product images:



Coomassie blue staining of purified ALOX15 protein (Cat# [TP306621]). The protein was produced from HEK293T cells transfected with ALOX15 cDNA clone (Cat# [RC206621]) using MegaTran 2.0 (Cat# [TT210002]).