

Product datasheet for **TP306630L**

FAM76A (NM_152660) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 76, member A (FAM76A), transcript variant 3, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206630 protein sequence Red =Cloning site Green =Tags(s) MAALYACTKCHQRFPEALSQGQQLCKECRIHPVWKCTYCRTEYQQESKTNICKKCAQNVQLYGTPKP CQYCNIIAAFIGNKCQRCTNSEKKYGPPYSCEQCKQQCAFDRKDDRKKVDGKLLCWLCTLSYKRVLQKTK EQRKHLSSSRAGHQEKEQYSRLSGGGHYNSQKTLSTSSIQNEIPKKKSKFESITNGDSFSPDLALDSP GTDHFVIAQLKEEVATLKKMLHQKQDMILEKEKKITELKADFQYQESQMRKMNQMEKTHKEVTEQLQ A KNRELLKQAAALSKSKKSEKSGAITSP TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	34.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_689873</u>



[View online »](#)

Locus ID: 199870
 UniProt ID: [Q8TAV0](#)
 RefSeq Size: 3360
 Cytogenetics: 1p35.3
 RefSeq ORF: 921

Product images:


Coomassie blue staining of purified FAM76A protein (Cat# [TP306630]). The protein was produced from HEK293T cells transfected with FAM76A cDNA clone (Cat# [RC206630]) using MegaTran 2.0 (Cat# [TT210002]).