

## Product datasheet for TP306650L

### CCDC89 (NM\_152723) Human Recombinant Protein

#### Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human coiled-coil domain containing 89 (CCDC89), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206650 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	<p>MRAPMLQKQQAPRMDTPPPEERLEKQNEKLNQEEETEFKELDGLREALANLRLGSEEERSEKAMLRSLR EEQSQLICILKRRSDEALERCQILELLNAELEEKMMQEAELKKAQGEYSRKLEERFMTLANHELMRLFK DEYKSENIKLRENEKLRLENSSLFSQALKDDEAKVLQLTVRCEALTGELETLKERCAQDACQAQAREKE LLELSQQAQACTHTKETEQLRSQLQTLKQQHQQAVEQIAKAEETHSSLSQELQARLQTVTREKEELLQSLI ERGKVLQNKQAEICQLEEKLEIANEDRKHALERFEQEAVALVDSNLRVRELQRKVDGIQKAYDELRLQSEA FKKHSLDLLSKERELNGKLRHLSP</p> <p><b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b></p>
Tag:	C-Myc/DDK
Predicted MW:	43.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_689936</a>
Locus ID:	220388



[View online »](#)

UniProt ID: [Q8N998](#)  
RefSeq Size: 2428  
Cytogenetics: 11q14.1  
RefSeq ORF: 1122

### Product images:



Coomassie blue staining of purified CCDC89 protein (Cat# [TP306650]). The protein was produced from HEK293T cells transfected with CCDC89 cDNA clone (Cat# [RC206650]) using MegaTran 2.0 (Cat# [TT210002]).