

Product datasheet for TP306657L

TMEM252 (NM_153237) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 9 open reading frame 71 (C9orf71), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206657 representing NM_153237 Red =Cloning site Green =Tags(s) MQNRTGLILCALALLMGFLMVCLGAFFISWGSIFDCQGSLIAAYLLLPLGFVILLSGIFWSNYRQVTESK GVLRHMLRQHLAGALPVATVDRPDFYPPAYEESLEVEKQSCPAEREASGIPPLYTETGLEFQDGNDSH PEAPPSYRESIAGLVVTAISEDAQRRGQEC TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	18.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_694969
Locus ID:	169693
UniProt ID:	Q8N6L7
RefSeq Size:	1284



[View online »](#)

Cytogenetics: 9q21.11
RefSeq ORF: 510
Synonyms: C9orf71
Protein Families: Transmembrane

Product images:

Coomassie blue staining of purified TMEM252 protein (Cat# [TP306657]). The protein was produced from HEK293T cells transfected with TMEM252 cDNA clone (Cat# [RC206657]) using MegaTran 2.0 (Cat# [TT210002]).