

Product datasheet for **TP306658L**

C12orf26 (METTL25) (NM_032230) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 12 open reading frame 26 (C12orf26), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206658 protein sequence Red =Cloning site Green =Tags(s)
	MAASCPLPVTPDLPTLRACLQGLLQFLRDALSISNAHTVDFYTESVWEELVDLPPEVLAAALRKSASETE ALPSETRPLVEAEWEAGMTDFPKIFCETSQKLVSVEAFALAAKYYSVQNLGICTPFEQLLVALRGNQNQR IGENQKAVEFMNMKKSHEVQAMSELISSADYYGIKQVIDLGSGKGYLSSFLSKYGLKVYIDSSNTNT HGAEERNRKLKKHWKLCHAQSRDLVNGLALKMAKERKVKNKVKNKADTEEVFNNSPTNQEKMPSTAIL PD FSGSVISNIRNQMETLHSQPHQEENLCFENSFSLINLLPINAVEPTSSQQIPNRETSEANKERRKMTSKS SESNIYSPLTSFITADSELHDIKDLEDCLMVGLHTCGDLAPNTRLRIFTSNSEIKGVCSVGCYHLLSEE FENQHKERTQEKWGFPMCHYLKEERWCCGRNARMSACLALERVAAGQGLPTESLFYRAVLQDIKDCYGI TKCDRHVGKIYSKCSSFLDYVRRSLKKLGLDESKLPEKIIMNYYEKYKPRMNELEAFNMLKVVLPACIET LILLDRLCYLKEQEDIAWSALVKLFDPVKSPRCYAVIALKKQQ TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	68 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_115606</u>
Locus ID:	84190
UniProt ID:	<u>Q8N6Q8</u>
RefSeq Size:	2105
Cytogenetics:	12q21.31
RefSeq ORF:	1809
Synonyms:	C12orf26
Summary:	Putative methyltransferase.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified METTL25 protein (Cat# [TP306658]). The protein was produced from HEK293T cells transfected with METTL25 cDNA clone (Cat# [RC206658]) using MegaTran 2.0 (Cat# [TT210002]).