

Product datasheet for **TP306664L**

MGC34774 (NM_203308) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human ribosomal protein L13A-like (MGC34774), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206664 representing NM_203308 Red =Cloning site Green =Tags(s) MNSVSVMMMLTSKPLHILGRFITKPLLFTNSMENFELTPSRHGRVHHQSCHEAPRTYKVAEDSRKTEPG AQWPRPSCGPPGSHLGQVGSAGKETGGRAPRGHQHFWQFLEKQIKVSGPLQADEHQPFPRPARPAIP GP PAHLLVALQNQVGPGCPGPPYGVLWDLTAL TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	19
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_976053</u>
Locus ID:	399670
RefSeq Size:	1848




[View online »](#)

Cytogenetics: 7q21.11

RefSeq ORF: 510

Product images:



Coomassie blue staining of purified MGC34774 protein (Cat# [TP306664]). The protein was produced from HEK293T cells transfected with MGC34774 cDNA clone (Cat# [RC206664]) using MegaTran 2.0 (Cat# [TT210002]).