

Product datasheet for TP306672

C11orf42 (NM_173525) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 11 open reading frame 42 (C11orf42), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206672 protein sequence Red =Cloning site Green =Tags(s)
	MLVGTPLLTLDEADATWTLIKDKVIEEHFGPNAVAVPFLSDAACDYLGLVVKQSRPAHTRLALPGRQG RRALKPVGPLPSLLEQAGSEGAFAHCTREYSPNGRAERAYEETRMLDGGQPCIRLHMGLRKKVAFLLLP PGQVSLQQTLPWLRSTHSIYVIYQVFSCSWLQLGLTSTAREPQLLRLLRSLPVAFSLKFSLQSKGVLGP QKPLTKDPLPHGANWVRPNLSIMPPLAPTSASADTTEAADVPPPVPAPPTPPPQEGPEDKPTRFSYKGRN PFWRGPQILSENWLFSPRSPPPGAQGGGPRDPDGHSMSLPLLQGLSSEFDSDD
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	36.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_775796</u>
Locus ID:	160298



[View online »](#)

UniProt ID: [Q8N5U0](#)
RefSeq Size: 1149
Cytogenetics: 11p15.4
RefSeq ORF: 999

Product images:



Coomassie blue staining of purified C11orf42 protein (Cat# TP306672). The protein was produced from HEK293T cells transfected with C11orf42 cDNA clone (Cat# [RC206672]) using MegaTran 2.0 (Cat# [TT210002]).