

Product datasheet for **TP306684**

C1orf87 (NM_152377) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human chromosome 1 open reading frame 87 (C1orf87), 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC206684 protein sequence
Red=Cloning site **Green**=Tags(s)

MSSAWKTPRGSDAMPEIMVKIIGSKHFQYLVEKPKIKENDSLKTETQTMHQKPMTDNARQMSRDTPVPIN
FTDQQTDDNPDDVKEKKHPENNQKSENNQKLLTGANSSRFLDGNIPSQANVHCSSVPTGDQSLSYVHGIP
RRKL RDWSLEQMVRGSSDQPEDIGQSPSGTTNEDAFLLALVRRELKSRPLSSNLEKLQKELKILDPISS
GFLLSQLSRLFLKHEVPLQLPTVKILCQRFSGKSPGEMVNYEKLLWFLNSAASDYPQQNKAADLRKTE
SHGTHSQSTPPQHSSSQPEVNRSLLEILKMA LR TTNGRLNIDNLNLSFRKEDRSFSGCLPLPKVRAICGK
HGLYLTLSLLETLLNHQDLGYQNEIKWQNFVEMLTRASSDLLSDLP TGKNEKKAPAPMEPEVPEMSQSK
TEHMKTP EEELQPESSPAETSACKDPLKPLKIRPVSQPFVNP AVKNKAAECETWIDRFKLENALYLCDL
SNTGVLEKERARRLIHNYNLIYNLSLSPQKIDQALRRFRSGENMLLEPALRYLKE L

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 61.9 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.


Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.



[View online »](#)

RefSeq: NP_689590
 Locus ID: 127795
 UniProt ID: Q8N0U7
 RefSeq Size: 2126
 Cytogenetics: 1p32.1
 RefSeq ORF: 1638
 Synonyms: CREF

Product images:


Coomassie blue staining of purified C1orf87 protein (Cat# TP306684). The protein was produced from HEK293T cells transfected with C1orf87 cDNA clone (Cat# [RC206684]) using MegaTran 2.0 (Cat# [TT210002]).