

Product datasheet for **TP306693L**

C18orf56 (TYMSOS) (NM_001012716) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 18 open reading frame 56 (C18orf56), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206693 protein sequence Red =Cloning site Green =Tags(s)
	MTPASGATASLGRLRARPRSRWDAAYLPAVAAVCVARASHVPNGTLRFGVCKARRTMRPLPGRIEVRTKR GPQRPAAPERSQPRLPPSRHPSRRGPRRHLSGCSAPACRIPTGCRCPCGRPS
	TR TRPLE QKLISEEDLA ANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_001012734
Locus ID:	494514
UniProt ID:	Q8TAI1
RefSeq Size:	860
Cytogenetics:	18p11.32




[View online »](#)

RefSeq ORF: 369

Synonyms: C18orf56

Product images:



Coomassie blue staining of purified TYMSOS protein (Cat# [TP306693]). The protein was produced from HEK293T cells transfected with TYMSOS cDNA clone (Cat# [RC206693]) using MegaTran 2.0 (Cat# [TT210002]).