

Product datasheet for TP306700


CCDC172 (NM_198515) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 10 open reading frame 96 (C10orf96)
Species:	Human
Expression Host:	HEK293T
Tag:	C-Myc/DDK
Predicted MW:	30.9 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_940917
Locus ID:	374355
RefSeq Size:	1703
Cytogenetics:	10q25.3
RefSeq ORF:	774
Synonyms:	C10orf96



[View online »](#)

Product images:

Coomassie blue staining of purified CCDC172 protein (Cat# TP306700). The protein was produced from HEK293T cells transfected with CCDC172 cDNA clone (Cat# [RC206700]) using MegaTran 2.0 (Cat# [TT210002]).