

Product datasheet for **TP306729M**


EAN57 (TEX33) (NM_178552) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 22 open reading frame 33 (C22orf33), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206729 protein sequence Red =Cloning site Green =Tags(s)
	MDPQSRSLKNAGSRSSSRENRATSGEGAQPCQGTDDGPSLGAQDQRSTPTNQKGSIIIPNNIRHKFGSNV V DQLVSEEQAQKAIDEVFEGQKRASSWPSRTQNPVEISSVFSDDYDLGYNMRSNLFRGAAEETKSLMKASY TPEIEKSVRDLEHWHGRKTDDLGRWHQKNAMNLLQKALEEKYGENSKSKSSKY
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	21.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_848647</u>
Locus ID:	339669
UniProt ID:	<u>O43247</u>


[View online »](#)

RefSeq Size:	1030
Cytogenetics:	22q12.3
RefSeq ORF:	585
Synonyms:	C22orf33; cE81G9.2; EAN57

Product images:


Coomassie blue staining of purified TEX33 protein (Cat# [TP306729]). The protein was produced from HEK293T cells transfected with TEX33 cDNA clone (Cat# [RC206729]) using MegaTran 2.0 (Cat# [TT210002]).