

Product datasheet for **TP306747**

MS4A15 (NM_152717) Human Recombinant Protein

Product data:


| | |
|---------------------------------------|---|
| Product Type: | Recombinant Proteins |
| Description: | Recombinant protein of human membrane-spanning 4-domains, subfamily A, member 15 (MS4A15), transcript variant 2, 20 µg |
| Species: | Human |
| Expression Host: | HEK293T |
| Expression cDNA Clone or AA Sequence: | >RC206747 protein sequence Red =Cloning site Green =Tags(s) |
| | MVRRGHVGIFFIEGGVPFWGGACFIISGSLSVAAEKNHTSCLVRSSLGTNLSVMAAFAGTAILLMDFGV TNRDVDRLGYLAVLTIFTVLEFFTAVIAMHFGCQAIHAQASAPVIFLPNAFSADFNIPSPAASAPPAYDNV AYAQGVV |
| | TR TRPLEQKLISEEDLAANDILDYKDDDDKV |
| Tag: | C-Myc/DDK |
| Predicted MW: | 15.3 kDa |
| Concentration: | >0.05 µg/µL as determined by microplate BCA method |
| Purity: | > 80% as determined by SDS-PAGE and Coomassie blue staining |
| Buffer: | 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol |
| Preparation: | Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps. |
| Note: | For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process. |
| Storage: | Store at -80°C. |
| Stability: | Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles. |
| RefSeq: | <u>NP_689930</u> |
| Locus ID: | 219995 |
| UniProt ID: | <u>Q8N5U1</u> |



[View online »](#)

| | |
|-------------------|--|
| RefSeq Size: | 1693 |
| Cytogenetics: | 11q12.2 |
| RefSeq ORF: | 441 |
| Summary: | May be involved in signal transduction as a component of a multimeric receptor complex. [UniProtKB/Swiss-Prot Function] |
| Protein Families: | Druggable Genome, Transmembrane |

Product images:



Coomassie blue staining of purified MS4A15 protein (Cat# TP306747). The protein was produced from HEK293T cells transfected with MS4A15 cDNA clone (Cat# [RC206747]) using MegaTran 2.0 (Cat# [TT210002]).