

Product datasheet for TP306785

PRR15 (NM_175887) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human proline rich 15 (PRR15), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC206785 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)
	MADSGDAGSSGPWWKSLTNSRKKSKEAAVGVPPPAQPAPGEPTPPAPPSPDWTSSSRENQHPNLLGGAGE PPKPDKLYGDKSGSSRRNLKISRSGRFKEKRKVRATLLPEAGRSPEEAGFGDPHEDKQ
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	13.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_787083
Locus ID:	222171
UniProt ID:	Q8IV56
RefSeq Size:	1717
Cytogenetics:	7p14.3


[View online »](#)

RefSeq ORF: 387

Summary: May have a role in proliferation and/or differentiation.[UniProtKB/Swiss-Prot Function]

Product images:



Coomassie blue staining of purified PRR15 protein (Cat# TP306785). The protein was produced from HEK293T cells transfected with PRR15 cDNA clone (Cat# [RC206785]) using MegaTran 2.0 (Cat# [TT210002]).