

## Product datasheet for **TP306790M**

### C10orf111 (NM\_153244) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 10 open reading frame 111 (C10orf111), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206790 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MESLQTPQHRENQDKREKEYGVKHMPMGNNAGNLEPEKRVAVRVALSSATAAQNIPSSVHCGCSKQWRLR LPSESLSRQVMKRPNNILKLRNLDLLIYPWPELRRRQVASDLMSLLLLPAFSGLTWAPFLFLFTYLPP FLNLLTVGFVSYFLV
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	17.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_694976</a>
Locus ID:	221060
RefSeq Size:	1785
Cytogenetics:	10p13




[View online »](#)

RefSeq ORF: 465

Synonyms: bA455B2.4

Protein Families: Transmembrane

### Product images:



Coomassie blue staining of purified C10orf111 protein (Cat# [TP306790]). The protein was produced from HEK293T cells transfected with C10orf111 cDNA clone (Cat# [RC206790]) using MegaTran 2.0 (Cat# [TT210002]).