

Product datasheet for **TP306815**

MRFAP1 (NM_033296) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Mof4 family associated protein 1 (MRFAP1), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206815 protein sequence Red =Cloning site Green =Tags(s)
	MRPLDIVELAEPEEVEVLEPEEDFEQFLLPVINEMREDIASLTREHGRAYLRNRSKLVEMDNMLIQIKTQ VEASEESALNHLQNPGDAAEGRAAKRCEKAEEKAKEIAKMAEMLVELVRRIEKSESS
	TR TRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	14.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_150638
Locus ID:	93621
UniProt ID:	Q9Y605 , B3KQA0
RefSeq Size:	2184
Cytogenetics:	4p16.1




[View online »](#)

RefSeq ORF: 381

Synonyms: PAM14; PGR1

Summary: This gene encodes an intracellular protein that interacts with members of the MORF4/MRG (mortality factor on chromosome 4/MORF4 related gene) family and the tumor suppressor Rb (retinoblastoma protein.) The protein may play a role in senescence, cell growth and immortalization. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Jan 2013]

Product images:



Coomassie blue staining of purified MRFAP1 protein (Cat# TP306815). The protein was produced from HEK293T cells transfected with MRFAP1 cDNA clone (Cat# [RC206815]) using MegaTran 2.0 (Cat# [TT210002]).