

Product datasheet for **TP306825L**

SPATA20 (NM_022827) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human spermatogenesis associated 20 (SPATA20), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA	>RC206825 protein sequence
Clone or AA Sequence:	Red=Cloning site Green=Tags(s)

MLGARAWLGRVLLLPRAGAGLAASRRCPGVWPRTWPHRSPSRGSSSRDKDRSATVSSSVMPAGGKGSHP
SSTPQRVPNRLIHEKSPYLLQHAYNPVDWYPWGQEAFDKARKENKPIFLSVGYSTCHWCHMMEESFQNE
EIGRLLSEDFVSVKVDREERPDVDKVMYTFVQATSSGGGWPMNVWLTPNLQPFVGGTYFPPEDGLTRVGF
RTVLLRIREQWKQNKNTLLENSQRVTTALLARSEISVGDRQLPPSAATVNNRCFQQLDEGYDEEYGGFAE
APKFPTPVILSFLFSYWLSHRLTQDGSRAQQMALHTLKMANGGIRDHVGGFHRYSTDRQWHVPHFEKM
LYDQAQLAVAYSQAFQLSGDEFYSVAKGILQYVARSLSHRSGGFYSAEDADSPPERGQRPKEGAYVWWT
VKEVQQLPEPVLGATEPLTSGQLLMKHGYLTEAGNISPSQDPKGELQGQNVLTVRYSELTAAARFGLDV
EAVRLLNSGLEKLFQARKHRPKPHLDSKMLAAWNGLMVSGYAVTGAVLGQDRLINATNGAKFLKRHMF
DVASGRMLMRTCYTGPGGTVEHSNPPCWGFLEDYAFVVRGLLDLYEASQESAWLEWALRLQDTQDKLFWDS
QGGGYFCSEAELGAGLPLRLKDDQDGAEPSANSVSAHNLLRLHGFTGHKDWMDKCVCLLTAFSERMRRVP
VALPEMVRALSAQQTLKQVICGDRQAKDTKALVQCVHSVYIPNKVLILADGDPSSFLSRQLPFLSTLR
RLEDQATAYVCENQACSVIPITDPCELRLLHP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	89.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_073738
Locus ID:	64847
UniProt ID:	Q8TB22
RefSeq Size:	2752
Cytogenetics:	17q21.33
RefSeq ORF:	2406
Synonyms:	HEL-S-98; SSP411; Tisp78
Summary:	May play a role in fertility regulation.[UniProtKB/Swiss-Prot Function]
Protein Families:	Druggable Genome

Product images:



Coomassie blue staining of purified SPATA20 protein (Cat# [TP306825]). The protein was produced from HEK293T cells transfected with SPATA20 cDNA clone (Cat# [RC206825]) using MegaTran 2.0 (Cat# [TT210002]).