

Product datasheet for **TP306829M**

TBX3 (NM_016569) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human T-box 3 (TBX3), transcript variant 2, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206829 protein sequence Red =Cloning site Green =Tags(s)

MSLSMRDPVIPGTSMAYHPFLPHRAPDFAMSAVLGHQPPFFPALTLPNGAAALSLPGALAKPIMDQLVG
AAETGIPFSSLGPQAHLRPLKTMPEEEEEVEDDPKVHLEAKELWDQFHKRGTEMVITKSGRRMFPPFKVRC
SGLDKKAKYILLMDIIAADDCRYKFHNSRWMVAGKADPEMPKRMYPHPDSPATGEQWMSKVWTFHKLKL
NNISDKHGFTLAFPSDHATWQGNYSFGTQILNSMHKYQPRFHIVRANDILKLPYSTFRTYLFPETEFIA
VTAYQNDKITQLKIDNNPFAKGFRTDGNRREKRKQLTLQSMRVFDERHKKENGTSDSSSEQAAFNCFA
QASSPAASTVGTSNLKDLCPEGESDAEAEESKEEHGPEACDAAKISTTTSEPCRDKGSPAVKAHLFAAE
RPRDSGRDLKASPDSRHSPATISSSTRGLGAEERRSPVREGTAPAKVEEARALPGKEAFAPLTVQTDAAA
AHLAQGPLPGLGFAPLAGQOFFNGHPLFLHPSQFAMGGAFSSMAAAGMGPLLATVSGASTGVSGLDSTA
MASAAAQGLSGASAATLPFHLQQHVLASQGLAMSPFGSLFPYPYTYMAAAAAASSAAASSSVHRHPFLN
LNTMRPRLRYPYSIPVVPDGSLLTTALPSMAAAGPLDGKVAALAASPASVAVDGSELNSRSSTLS
SSMSLSPKLCAEKEAATSELQSIQRLVSGLEAKPDRSRASAP

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	79.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_057653
Locus ID:	6926
UniProt ID:	O15119 , A0A024RBL6
RefSeq Size:	4814
Cytogenetics:	12q24.21
RefSeq ORF:	2229
Synonyms:	TBX3-ISO; UMS; XHL
Summary:	This gene is a member of a phylogenetically conserved family of genes that share a common DNA-binding domain, the T-box. T-box genes encode transcription factors involved in the regulation of developmental processes. This protein is a transcriptional repressor and is thought to play a role in the anterior/posterior axis of the tetrapod forelimb. Mutations in this gene cause ulnar-mammary syndrome, affecting limb, apocrine gland, tooth, hair, and genital development. Alternative splicing of this gene results in three transcript variants encoding different isoforms; however, the full length nature of one variant has not been determined. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Transcription Factors

Product images:



Coomassie blue staining of purified TBX3 protein (Cat# [TP306829]). The protein was produced from HEK293T cells transfected with TBX3 cDNA clone (Cat# [RC206829]) using MegaTran 2.0 (Cat# [TT210002]).