

Product datasheet for **TP306846L**

RBFOX2 (NM_001031695) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins
Description: Purified recombinant protein of Homo sapiens RNA binding motif protein 9 (RBM9), transcript variant 1, 1 mg
Species: Human
Expression Host: HEK293T
Expression cDNA Clone or AA Sequence: >RC206846 protein sequence
Red=Cloning site **Green**=Tags(s)

MEKKKMVTQGNQEPTTTPDAMVQPFTTIPFPPPPQNGIPTEYGVPHQTQDYAGQTGEHNLTYGSTQAHGE
QSSNSPSTQNGSLTTEGGAQTDGQQSQTQSSENSESKSTPKRLHVSNIPIFRFRDPLRQMFGQFGKILDV
EIIENERGSKGFGFVTFENSADADRAREKLHGTVEGRKIEVNNATARVMTNKKMVTYPYANGWKLSPVVG
AVYGPELYAASSFQADVSLGNDAAVPLSGRGGINTYIPLISLPLVPGFPYPTAATTAFAFRGAHLRGRGR
TVYGAVRAVPPTAIPAYPGVVYQDGFYGDLYGGYAAARYAQPATATAATAAAAAAAAAAAYSDGYGRVYTAD
PYHALAPAASYGVGAVASLYRGGYSRFAPY

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK
Predicted MW: 40.2 kDa
Concentration: >0.05 µg/µL as determined by microplate BCA method
Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage: Store at -80°C.
Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq: [NP_001026865](#)




[View online »](#)

Locus ID: 23543
UniProt ID: [O43251](#)
RefSeq Size: 6997
Cytogenetics: 22q12.3
RefSeq ORF: 1140
Synonyms: dJ106l20.3; Fox-2; FOX2; fxx; HNRBP2; HRNBP2; RBM9; RTA

Summary: This gene is one of several human genes similar to the *C. elegans* gene Fox-1. This gene encodes an RNA binding protein that is thought to be a key regulator of alternative exon splicing in the nervous system and other cell types. The protein binds to a conserved UGCAUG element found downstream of many alternatively spliced exons and promotes inclusion of the alternative exon in mature transcripts. The protein also interacts with the estrogen receptor 1 transcription factor and regulates estrogen receptor 1 transcriptional activity. Multiple transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jul 2008]

Protein Families: Druggable Genome, Transcription Factors

Product images:



Coomassie blue staining of purified RBFOX2 protein (Cat# [TP306846]). The protein was produced from HEK293T cells transfected with RBFOX2 cDNA clone (Cat# [RC206846]) using MegaTran 2.0 (Cat# [TT210002]).