

Product datasheet for TP306852

ZNF239 (NM_005674) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human zinc finger protein 239 (ZNF239), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206852 protein sequence Red =Cloning site Green =Tags(s)
	<p>MASTITGSQDCIVNHRGEVDGPELDISPCQQWGEASSPISRNRDSVMTLQSGCFENIESETYLPLKVSS QIDTQDSSVKFCKNEPQDHQESRRLFVMEESTERKVIKGESCSENLQVKLVSDGQELASPLLNGEATCQN GQLKESLDPIDCNCKDIHGWSQVWSCSQRAHTEEKPCDHNNCGKILNTSPDGHPYEKIHTAEKQYEC QCGKNFSQSSELLHQRDHTTEKPYKCEQCGKGFTRSSLLIHQAVHTDEKPYKCDKCGKGFTRSSLLI HHAVHTGEKPYKCDKCGKGFSSKLHIIHQRVHTGEKPYECECGMSFSQRSNLHIIHQRVHTGERPYKCG ECGKGFSSNLHIIHRCIHTGEKPYQCYECGKGFSSDLRIHLRVHTGEKPYHCGKCGKGFSSKLLI HQRVHTGEKPYECSKCGKGFSSNLHIIHQRVHKKDPR</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	51.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_005665</u>



[View online »](#)

Locus ID:	8187
UniProt ID:	Q16600
RefSeq Size:	2406
Cytogenetics:	10q11.21
RefSeq ORF:	1374
Synonyms:	HOK-2; MOK2
Summary:	MOK2 proteins are DNA- and RNA-binding proteins that are mainly associated with nuclear RNP components, including the nucleoli and extranucleolar structures (Arranz et al., 1997 [PubMed 9121460]).[supplied by OMIM, Mar 2008]
Protein Families:	Transcription Factors