

Product datasheet for **TP306860**

HPD (NM_002150) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human 4-hydroxyphenylpyruvate dioxygenase (HPD)
Species:	Human
Expression Host:	HEK293T
Tag:	C-Myc/DDK
Predicted MW:	44.8 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_002141
Locus ID:	3242
RefSeq Size:	1440
Cytogenetics:	12q24.31
RefSeq ORF:	1179
Synonyms:	4-HPPD; 4HPPD; GLOD3; HPPDASE; PPD
Summary:	The protein encoded by this gene is an enzyme in the catabolic pathway of tyrosine. The encoded protein catalyzes the conversion of 4-hydroxyphenylpyruvate to homogentisate. Defects in this gene are a cause of tyrosinemia type 3 (TYRO3) and hawkinsinuria (HAWK). Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Jan 2010]
Protein Families:	Druggable Genome



[View online »](#)

Protein Pathways: Metabolic pathways, Phenylalanine metabolism, Tyrosine metabolism, Ubiquinone and other terpenoid-quinone biosynthesis

Product images:



Coomassie blue staining of purified HPD protein (Cat# TP306860). The protein was produced from HEK293T cells transfected with HPD cDNA clone (Cat# [RC206860]) using MegaTran 2.0 (Cat# [TT210002]).