

## Product datasheet for **TP306877**

### UBALD1 (NM\_145253) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human family with sequence similarity 100, member A (FAM100A), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206877 protein sequence <b>Red</b> =Cloning site <b>Green</b> =Tags(s)
	MSVNMDELKHQVMINQFVLTAGCAADQAKQLLQAAHWQFETALSAFFQETNIPYSHHHHQMMPANTPA TPPNFPDALTMFSRLKASESFHSGGSGSPMAATATSPPHFPHAATSSSAASSWPTAASPPGGPQHHQPQ PPLWTPTPPSPASDWPLAPQQATSEPRHPAMEAER
	<b>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</b>
Tag:	C-Myc/DDK
Predicted MW:	18.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_660296</a>
Locus ID:	124402
UniProt ID:	<a href="#">Q8TB05</a>
RefSeq Size:	1430



[View online »](#)

Cytogenetics: 16p13.3  
RefSeq ORF: 531  
Synonyms: FAM100A; PP11303

**Product images:**

Coomassie blue staining of purified UBALD1 protein (Cat# TP306877). The protein was produced from HEK293T cells transfected with UBALD1 cDNA clone (Cat# [RC206877]) using MegaTran 2.0 (Cat# [TT210002]).