

Product datasheet for **TP306919L**

PALM2AKAP2 (NM_001037293) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human paralemmin 2 (PALM2), transcript variant 2, 1 mg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC206919 protein sequence
Red=Cloning site Green=Tags(s)

MAEAE LHKERLQAI AEKRKRQTEIEGKRQQLDEQILLQHSKSKVLREK WLLQGIPAGTAE EEE EARRRQS
EEDEF RVKQLEDNIQRLEQEIQ TLESEESRISAKEQIILEKLKETEK SFKDFQKGF SSTDGAVYAMEINV
EKDKQTGETKILSTSTIGPEGVHQGVKVYDDG TKVVEVRS GGTVVENG VHKLSTKDVEELIQKAGQSS
LGGGHVSERTVIADGSL SHPK EHM LCKEAKLEMVHKS RDKHSSGNP GQQAQAPS AAGPEANLDQPVTMIF
MGYQNI EDEEETKKVLGYDETIK AELVLIDEDDEKSLREKT VTDVSTIDGNAAELVSGRPVSDTTEPSSP
EGKEESLATEPAPGTQKKKRCQCCVVM

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 42 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_001032370](#)


Locus ID: 445815



[View online »](#)

UniProt ID:	Q8IXS6
RefSeq Size:	9286
Cytogenetics:	9q31.3
RefSeq ORF:	1131
Synonyms:	AKAP-2; AKAP-KL; AKAP2; AKAPKL; MISP2; PALM2; PALM2-AKAP2; PRKA2
Summary:	This gene belongs to the paralemmin downstream gene (PDG) family defined in PMID:22855693. Paralemmin downstream genes may have evolved contiguously with the paralemmin genes and are associated with other paralemmin paralogs in humans and several other taxa. The gene encodes three distinct protein isoforms, the PALM2 isoform, the AKAP2 isoform and the PALM2-AKAP2 isoform. The biological significance of the PALM2-AKAP2 isoforms is yet unknown. Earlier, PALM2 and AKAP2 were annotated as separate genes and PALM2-AKAP2 was annotated as a readthrough gene. [provided by RefSeq, May 2019]

Product images:



Coomassie blue staining of purified PALM2-AKAP2 protein (Cat# [TP306919]). The protein was produced from HEK293T cells transfected with PALM2-AKAP2 cDNA clone (Cat# [RC206919]) using MegaTran 2.0 (Cat# [TT210002]).