

Product datasheet for TP306937L

C19orf34 (NM_152771) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 19 open reading frame 34 (C19orf34), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206937 representing NM_152771 <div> <div>Red</div>=Cloning site <div>Green</div>=Tags(s) </div> <p>MGRAGPRSTADTRPRATVITTRRPWQKPTSPRRLHRRRPRGQPAENAAEPSQVAISKGKSPDYANGG KAAMDSRASDAINKSKVDTSASNPSQRRPSRLRGKRSLALPPRLECLFPSSCGSRRTSRPGFPPLGSH LGCQLLPL</p> <p><div>TR</div><div>TRPLEQKLISEEDLAANDILDYKDDDDKV</div></p>
Tag:	C-Myc/DDK
Predicted MW:	15.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_689984</u>
Locus ID:	255193
RefSeq Size:	1281
Cytogenetics:	19p13.3


[View online »](#)

RefSeq ORF: 729

Synonyms: chromosome 19 open reading frame 34; MGC39696

Product images:



Coomassie blue staining of purified C19orf34 protein (Cat# [TP306937]). The protein was produced from HEK293T cells transfected with C19orf34 cDNA clone (Cat# [RC206937]) using MegaTran 2.0 (Cat# [TT210002]).