

Product datasheet for TP306957M

SEPTIN6 (NM_015129) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human septin 6 (SEPT6), transcript variant II, 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206957 protein sequence Red=Cloning site Green=Tags(s)
	<p>MAATDIARQVGEGCRTVPLAGHVGFDLSPDQLVNKSVSQGFCFNILCVGETGLGKSTLMDTLFNTKFEGE PATHTQPGVQLQSNTYDLQESNVRLKLTIVSTVGFQDQINKEDSYKPIVEFIDAQFEAYLQEELKIRRVL HTYHDSRIHVCLYFIAPTGHSLKSLDLVTMKKLDSKVNIIIAKADAISKSELTKFKIKITSELVSNV QIYQFPTDDESVAEINGTMNAHLPFAVIGSTEELKIGNKMMRARQYPWGTQVENEAHCDFVKLREMLIR VNMEDLREQTHTRHYELYRCKLEEMGFKDTPDPSKPFSLQETYEAKRNEFLGELQKKEEEMRQMFVQRV KEKEAELKEAEKELHEKFDRLKHLHQDEKKKLEDKKKSLDDEVNAFKQRKTA AELLQSQGSQAGGSQTLK RDKEKKNPWLCTE</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	49.5 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_055944</u>




View online »

Locus ID:	23157
UniProt ID:	Q14141
RefSeq Size:	3387
Cytogenetics:	Xq24
RefSeq ORF:	1302
Synonyms:	SEP2; SEPT2; SEPT6

Summary: This gene is a member of the septin family of GTPases. Members of this family are required for cytokinesis. One version of pediatric acute myeloid leukemia is the result of a reciprocal translocation between chromosomes 11 and X, with the breakpoint associated with the genes encoding the mixed-lineage leukemia and septin 2 proteins. This gene encodes four transcript variants encoding three distinct isoforms. An additional transcript variant has been identified, but its biological validity has not been determined. [provided by RefSeq, Jul 2008]

Product images:



Coomassie blue staining of purified SEPTIN6 protein (Cat# [TP306957]). The protein was produced from HEK293T cells transfected with SEPTIN6 cDNA clone (Cat# [RC206957]) using MegaTran 2.0 (Cat# [TT210002]).