

Product datasheet for TP306984

INPP5E (NM_019892) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human inositol polyphosphate-5-phosphatase, 72 kDa (INPP5E), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC206984 representing NM_019892 Red=Cloning site Green=Tags(s)

MPSKAENLRPSEPAPQPPEGRTLQGQLPGAPPAQRAGSPPDAPGSESPALACSTPATPSGEDPPARAAPI
APRPPARPRLERALSDDKGWRRRRFRGSQEDLEARNGTSPSRGVSQSEGGAPAHSCSPCLSTSLQEI
PKSRGVLSSERGSPSSGGNPLSGVASSPNLPHRDAAVAGSSPRLPSLLPPRPPPPALSLDIASDSLRTAN
KVDSLADYKLRAQPLLVAHSSLGPRPRSPPLACDDCSLRSKSSFSLLAPIRSKDVRSRSYLEGSLLA
SGALLGADELARYFPDRNVALFVATWNMQGQKELPPSLDEFLLPAEADYAQDLYVIGVQEGCSDRREWET
RLQETLGPYVLLSSAAHGVLYMSLFIRRDLIWFCSEVECSTVTRIVSQIKTKGALGISFTFFGTSLF
ITSHFTSGDGKVAERLLDYTRTVQALVLPNVPDTNPYRSSAADVTRFDEVFVWFGDFNRLSGGRTVVD
ALLCQGLVVDVPALLQHDQLIREMRKGSIFKGFQEPDIHFLPSYKFDIGKDYDSTSKQRTPSYTDRLVY
RSRHKGDICPVSYS SCPGIKTS DHRPVYGLFRVKVRPGRDNIPLAAGKFDRELYLLGIKRRISKEIQRQQ
ALQSQNSSTICSVS

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	70 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_063945
Locus ID:	56623
UniProt ID:	Q9NRR6
RefSeq Size:	3440
Cytogenetics:	9q34.3
RefSeq ORF:	1932
Synonyms:	CORS1; CPD4; JBTS1; MORMS; pharbin; PPI5PIV
Summary:	The protein encoded by this gene is an inositol 1,4,5-trisphosphate (InsP3) 5-phosphatase. InsP3 5-phosphatases hydrolyze Ins(1,4,5)P3, which mobilizes intracellular calcium and acts as a second messenger mediating cell responses to various stimulation. Studies of the mouse counterpart suggest that this protein may hydrolyze phosphatidylinositol 3,4,5-trisphosphate and phosphatidylinositol 3,5-bisphosphate on the cytoplasmic Golgi membrane and thereby regulate Golgi-vesicular trafficking. Mutations in this gene cause Joubert syndrome; a clinically and genetically heterogenous group of disorders characterized by midbrain-hindbrain malformation and various associated ciliopathies that include retinal dystrophy, nephronophthisis, liver fibrosis and polydactyly. Alternative splicing results in multiple transcript variants encoding different isoforms. [provided by RefSeq, Jan 2016]
Protein Families:	Druggable Genome
Protein Pathways:	Inositol phosphate metabolism, Metabolic pathways, Phosphatidylinositol signaling system

Product images:



Coomassie blue staining of purified INPP5E protein (Cat# TP306984). The protein was produced from HEK293T cells transfected with INPP5E cDNA clone (Cat# [RC206984]) using MegaTran 2.0 (Cat# [TT210002]).