

Product datasheet for **TP307057L**

C13orf26 (TEX26) (NM_152325) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 13 open reading frame 26 (C13orf26), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207057 protein sequence Red =Cloning site Green =Tags(s)
	MEQPGPRAPDPSLCHHNLQPTDDPNWDSYATTMRTAFTPKTGAVPALIRQNGIRRLGYTYSLSDPILNQ QYSDEYTWKSHSKEDLIKTETSRGIKSHKSHLNEDIFLWTLPHCQQTGTLKNCLPWKIPASMKEVNKALS NQFISLTKRDFVDRSKAQKIKKSSHLSLEWKLLPQPPDTEFRRNYQIPAKIPELQDFSFKYGCYSSLPV ASQGLVPSVLHSYLRNQEHTKKQTTYQSDYDKTYPDFLMLLSFTSSQVKEYLQSLSYKDRQIIDRFIRT HCDTNKKKK
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	33.4 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_689538</u>
Locus ID:	122046



[View online »](#)

UniProt ID:	<u>Q8N6G2</u>
RefSeq Size:	1113
Cytogenetics:	13q12.3
RefSeq ORF:	867
Synonyms:	C13orf26

Product images:


Coomassie blue staining of purified TEX26 protein (Cat# [TP307057]). The protein was produced from HEK293T cells transfected with TEX26 cDNA clone (Cat# [RC207057]) using MegaTran 2.0 (Cat# [TT210002]).