

Product datasheet for **TP307059M**

TEX9 (NM_198524) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human testis expressed 9 (TEX9), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207059 protein sequence Red =Cloning site Green =Tags(s)
	<p>MAGRSLCLTRSSVPGTPFPPPVQQPSTPGPDLLALEEEYKRLNAELQAKTADVVQQAKEIIRDRQEVRSR PVSTQMKSCDEDDYSLRGLLPSEGIVHLHSETPKTKNIDPVNVQVQNKLSANKGRKTNSVVKLYSDV QTADDVAIPEDFSDFLAKTISKIEGQLEEEGLPEYIDDI FSGVSN DIGTEAQIRFLKAKLHVMQEELDN VVCECNKKEDEIQNLKSQVKNFEEDFMRQQRTINMQSQVEKYKTLFEEANKKYDGLQQQLSSVERELE N KRRLLQQAASSQSATEVRLNRAL EEA EKYKLELSKLRQNNKD IANEEHKKIEVLKSENKKLEKQKGELMI GFKKQLKLIDVLKRQKM HIEAAKMLSFT EEEFMKALEWGN S</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	44.6 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_940926</u>



[View online »](#)

Locus ID: 374618
 UniProt ID: Q8N6V9
 RefSeq Size: 1559
 Cytogenetics: 15q21.3
 RefSeq ORF: 1173

Product images:


Coomassie blue staining of purified TEX9 protein (Cat# [TP307059]). The protein was produced from HEK293T cells transfected with TEX9 cDNA clone (Cat# [RC207059]) using MegaTran 2.0 (Cat# [TT210002]).