

Product datasheet for **TP307079**

CD99L2 (NM_031462) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human CD99 molecule-like 2 (CD99L2), transcript variant 1, 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone >RC207079 protein sequence

or AA Sequence: **Red**=Cloning site **Green**=Tags(s)

MVAWRS AFLVCLAFSLATLVQRGSGDFDDFNLEDAVKETSSVKQPWDHTTTTTTNRPGTTRAPAKPPGSG
LDLADALDDQDDGRRKPGIGGRERWNHVTTTTTKRPVTTTRAPANTLGNDFDLADALDDRNRDDGRRKPIA
GGGGFSDKDLEDIVGGGEYKPKDKGKGDGRYGSNDDPGSGMVAEPGTIAGVASALAMALIGAVSSYISYQQ
KKFCFSIQQLNADYVKGENLEAVVCEEPQVKYSTLHTQSAEPPPPPEPARI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 27.8 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_113650](#)

Locus ID: 83692

UniProt ID: [Q8TCZ2](#), [A0A024RC16](#)



[View online >](#)

RefSeq Size: 3736

Cytogenetics: Xq28


RefSeq ORF: 786

Synonyms: CD99B; MIC2L1

Summary: This gene encodes a cell-surface protein that is similar to CD99. A similar protein in mouse functions as an adhesion molecule during leukocyte extravasation. Alternate splicing results in multiple transcript variants. [provided by RefSeq, May 2010]

Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified CD99L2 protein (Cat# TP307079). The protein was produced from HEK293T cells transfected with CD99L2 cDNA clone (Cat# [RC207079]) using MegaTran 2.0 (Cat# [TT210002]).