

Product datasheet for **TP307108**

NDUFS6 (NM_004553) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human NADH dehydrogenase (ubiquinone) Fe-S protein 6, 13kDa (NADH-coenzyme Q reductase) (NDUFS6), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207108 protein sequence Red =Cloning site Green =Tags(s)
	MAAAMTFCRLLNRCGEAARSLPLGARCFGVRVSPTGEKVTHHTGQVYDDKDYRRIRFVGRQKEVNENFAID LIAEQPVSEVETRVIACDGGGGALGHPKVYINLDKETKTGTGCGYCGLQFRQHHH
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	10.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_004544
Locus ID:	4726
UniProt ID:	O75380 , Q6IBC4
RefSeq Size:	554



[View online »](#)

Cytogenetics: 5p15.33


RefSeq ORF: 372

Synonyms: CI-13kA; CI-13kD-A; CI13KDA; MC1DN9

Summary: This gene encodes a subunit of the NADH:ubiquinone oxidoreductase (complex I), which is the first enzyme complex in the electron transport chain of mitochondria. This complex functions in the transfer of electrons from NADH to the respiratory chain. The subunit encoded by this gene is one of seven subunits in the iron-sulfur protein fraction. Mutations in this gene cause mitochondrial complex I deficiency, a disease that causes a wide variety of clinical disorders, including neonatal disease and adult-onset neurodegenerative disorders.[provided by RefSeq, Oct 2009]

Protein Pathways: Alzheimer's disease, Huntington's disease, Metabolic pathways, Oxidative phosphorylation, Parkinson's disease

Product images:



Coomassie blue staining of purified NDUFS6 protein (Cat# TP307108). The protein was produced from HEK293T cells transfected with NDUFS6 cDNA clone (Cat# [RC207108]) using MegaTran 2.0 (Cat# [TT210002]).