

Product datasheet for TP307230

CACNB1 (NM_000723) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human calcium channel, voltage-dependent, beta 1 subunit (CACNB1), transcript variant 1
Species:	Human
Expression Host:	HEK293T
Tag:	C-Myc/DDK
Predicted MW:	65.5 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<u>NP_000714</u>
Locus ID:	782
UniProt ID:	<u>Q02641</u>
RefSeq Size:	3687
Cytogenetics:	17q12
RefSeq ORF:	1794
Synonyms:	CAB1; CACNLB1; CCHLB1



[View online »](#)

Summary:


The protein encoded by this gene belongs to the calcium channel beta subunit family. It plays an important role in the calcium channel by modulating G protein inhibition, increasing peak calcium current, controlling the alpha-1 subunit membrane targeting and shifting the voltage dependence of activation and inactivation. Alternative splicing occurs at this locus and three transcript variants encoding three distinct isoforms have been identified. [provided by RefSeq, Jul 2008]

Protein Families:

Druggable Genome, Ion Channels: Other

Protein Pathways:

Arrhythmogenic right ventricular cardiomyopathy (ARVC), Cardiac muscle contraction, Dilated cardiomyopathy, Hypertrophic cardiomyopathy (HCM), MAPK signaling pathway

Product images:

Coomassie blue staining of purified CACNB1 protein (Cat# TP307230). The protein was produced from HEK293T cells transfected with CACNB1 cDNA clone (Cat# [RC207230]) using MegaTran 2.0 (Cat# [TT210002]).