

Product datasheet for **TP307285**

PLEKHA1 (NM_021622) Human Recombinant Protein

Product data:

Product Type: Recombinant Proteins

Description: Recombinant protein of human pleckstrin homology domain containing, family A (phosphoinositide binding specific) member 1 (PLEKHA1), transcript variant 1, 20 µg

Species: Human

Expression Host: HEK293T

Expression cDNA Clone or AA Sequence: >RC207285 protein sequence
Red=Cloning site **Green**=Tags(s)

MPYVDRQNRICGFLDIEENENSGKFLRRYFILDTREDSFVWYMDNPQNLPSGSSRVGAIKLTYSKVS
TKLRPKAEFCFVMNAGMRKYFLQANDQQDLVEWVNVLNKAIKITVPKQSDSQPNSDNL
SRHGECGKKQVS
YRTDIVGGVPIITPTQKEEVNECGESIDRNNLKRSQSHLPYFTPKPPQDSAVIKAGYCVKQ
GAVMKNWKR
RYFQLDENTIGYFKSELEKEPLRVIPLKEVHKVQECKQSDIMMRDNLFEIVTTSRTFYVQ
ADSPEEMHSW
IKAVSGAIVAQRGPGRSASSEHPPGPSESKHAFRPTNAATATSHSTASRSNSLSTFTMEK
RGFYESLAK
VKPGNFVKVQTVSPREPASKVTEQALLRPQSKNGPQEKDCDLVDLDDASLPVSDV

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag: C-Myc/DDK

Predicted MW: 45.4 kDa

Concentration: >0.05 µg/µL as determined by microplate BCA method

Purity: > 80% as determined by SDS-PAGE and Coomassie blue staining

Buffer: 25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol

Preparation: Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.

Note: For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.

Storage: Store at -80°C.

Stability: Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.

RefSeq: [NP_067635](#)




[View online >](#)

Locus ID: 59338
UniProt ID: [Q9HB21](#)
RefSeq Size: 3732
Cytogenetics: 10q26.13
RefSeq ORF: 1212
Synonyms: TAPP1

Summary: This gene encodes a pleckstrin homology domain-containing adapter protein. The encoded protein is localized to the plasma membrane where it specifically binds phosphatidylinositol 3,4-bisphosphate. This protein may be involved in the formation of signaling complexes in the plasma membrane. Polymorphisms in this gene are associated with age-related macular degeneration. Alternate splicing results in multiple transcript variants. A pseudogene of this gene is found on chromosome 5.[provided by RefSeq, Sep 2010]

Product images:



Coomassie blue staining of purified PLEKHA1 protein (Cat# TP307285). The protein was produced from HEK293T cells transfected with PLEKHA1 cDNA clone (Cat# [RC207285]) using MegaTran 2.0 (Cat# [TT210002]).