

## Product datasheet for TP307300L

### C7orf61 (NM\_001004323) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human hypothetical LOC402573 (LOC402573), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207300 protein sequence <span style="color: red;">Red</span> =Cloning site <span style="color: green;">Green</span> =Tags(s)  MVVVMKFFRWVRRRAWQRIISWVFFWRQKIKPTISGHPDSKKHSLKKMEKTLQVVETLRLVELPKAKPKL GESPELADPCVLAKTTEETEVELGQQGQSLLQLPRTAVKSVSTLMVSALQSGWQMCWKSSVSSASVSSQ VRTQSPLKTPEAELLWEVYLVLVAVRKHLRRLYRRQERHRRHHVRCHAAPRPNPAQSLKLDASPL  <span style="color: red;">TR</span> <span style="color: green;">TRPLEQKLISEEDLAANDILDYKDDDDKV</span>
Tag:	C-Myc/DDK
Predicted MW:	23.7 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001004323</a>
Locus ID:	402573
UniProt ID:	<a href="#">Q8IZ16</a>
RefSeq Size:	1020


[View online »](#)

Cytogenetics: 7q22.1

RefSeq ORF: 618

### Product images:



Coomassie blue staining of purified C7orf61 protein (Cat# [TP307300]). The protein was produced from HEK293T cells transfected with C7orf61 cDNA clone (Cat# [RC207300]) using MegaTran 2.0 (Cat# [TT210002]).