

Product datasheet for **TP307305L**

SMYD4 (NM_052928) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human SET and MYND domain containing 4 (SMYD4), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207305 protein sequence Red =Cloning site Green =Tags(s)
	<p>MDLPVDEWKSYPQLQKQWASLPTSVQVTISTAETLRDIFLHSSLLQPEDELFLKRLSKGYLVGKDSAPLF YREEGNKKFQEKDYTGAAVLYSKGVSHSRPNTEDMSLCHANRSAALFHLGQYETCLKDINIAQTHGYPER LQPKIMLRKAECLVALGRLQEASQTISDLERNFTATPALADVLPQTLQRNLHRLKMKMQEKDSLTFSPA ALAKTLEDAALREENEQLSNASSSIGLCVDPLKGRCLVATKDILPGELLVQEDAFVSVLNPGLPPP HHG LDSKWDTRVTNGDLYCHRCLKHTLATVPCDGC SYAKYCSQECLQQVWELYHRTECPLGGLLLTLGVFCHI ALRLTLLVGFEDVRKIITKCLKDKISNKDICLPESNNQVKTLNYGLGESEKNGNIVETPIPGCDINGKYEN NYNAVFNLLPHTENHSPEHKFLCALCVSALCRQLEAASLQAIPTERIVNSSQLKAAVTPELCPDVTI WGV AMLRHMLLQQCNAQAMTTIQHTGPKGSIVTDSRQVRLATGIFPVISLLNHSCSPNTSVSFISTVATIRAS QRIRKGQEILHCYGP HKSRMGVAERQQKLR SQYFFDCAC PACQTEAHRMAAGPRWEAFCCNSCGAPMQGD DVLRCGSRSCAESAVSRDHLVSRDLQDLQQQVRVAQKLLRDGELERAVQRLSGCQRDAESFLWAEHAVGE IADGLARACAALGDWQKSATHLQRSLCVVEVRHGPSSVEMGHELKLAQIFFNGFAVPEALSTIQKAEV LSLHCGPWDDEIQELQKMKSCLLDLPPTVGPAL</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	89 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.



[View online »](#)

Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_443160
Locus ID:	114826
UniProt ID:	Q8IYR2
RefSeq Size:	4349
Cytogenetics:	17p13.3
RefSeq ORF:	2412
Synonyms:	ZMYND21

Product images:



Coomassie blue staining of purified SMYD4 protein (Cat# [TP307305]). The protein was produced from HEK293T cells transfected with SMYD4 cDNA clone (Cat# [RC207305]) using MegaTran 2.0 (Cat# [TT210002]).