

## Product datasheet for TP307309L

### C2orf52 (NM\_173513) Human Recombinant Protein


#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human chromosome 2 open reading frame 52 (C2orf52), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207309 representing NM_173513 <span style="color: red;">Red</span> =Cloning site <span style="color: green;">Green</span> =Tags(s)  MSEAKDNGSRDEVLPVPHKNCRKNTTVPGKKGEEKSLAPVFAEKLISPSRRGAKLKDRESHQENEDRNSEL DQDEEDKESFCRGFPMMSGCELETSCCVCHSTALGERFC  <span style="color: red;">TR</span> TRPLE <span style="color: green;">QKLISEEDLA</span> AND <span style="color: green;">ILDYKDDDDK</span> V
Tag:	C-Myc/DDK
Predicted MW:	11.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_775784</a>
Locus ID:	151477
RefSeq Size:	1299
Cytogenetics:	2q37.1
RefSeq ORF:	324


[View online »](#)

Synonyms: C2orf52

### Product images:



Coomassie blue staining of purified C2orf52 protein (Cat# [TP307309]). The protein was produced from HEK293T cells transfected with C2orf52 cDNA clone (Cat# [RC207309]) using MegaTran 2.0 (Cat# [TT210002]).