

Product datasheet for TP307337

CRY2 (NM_021117) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cryptochrome 2 (photolyase-like) (CRY2), transcript variant 1, 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207337 protein sequence Red =Cloning site Green =Tags(s)

MAATVATAAAVAPAPAPGTDSASSVHWFRKGLRLHDNPALLAAVARGARCVRCVYILDPWFAASSSSVGINR
WRFLQSLLEDLDTSLRKLNSRFLVVRGQPADVFPRLFKEWGVTRLTFFEYDSEPFGKERDAAIMKMAKEAG
VEVVTENSHTLYDLDRIIELNGQKPLTYKRFQAIISRMELPKKPVGLVTSQQMESCRAEIQENHDETYG
VPSLEELGFPTTEGLGPAVWQGGETEALARLDKHLERKAWVANYERPRMNANSLLASPTGLSPYLRFGCLS
CRLFYYRLWDLYKKVKRNSTPPLSLFGQLLWREFFYTAATNNPRFDRMEGNPICIQIPWDRNPEALAKWA
EGKTGFPWIDAIMTQLRQEGWIHHLARHAVACFLTRGDLWVSWESGVRVFEDELLDADFSVNAGSWMWLS
CSAFFQQFFHCYCPVGFGRRTDPSGDYIRRYLPKLFKAFPSRYIYEPWNPESIQKAAKCIIGVDYPRPIV
NHAETSRNLNIERMKQIYQQLSRYRGLCLLASVPSCVEDLSHPVAEPSSSQAGSMSSAGPRPLPSGPASPK
RKLEAAEPPGEELSKRARVAELPTPELPSKDA

TRTRPLEQKLISEEDLAANDILDYKDDDDKV


Tag:	C-Myc/DDK
Predicted MW:	66.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_066940
Locus ID:	1408
UniProt ID:	Q49AN0 , A0A0D2X7Z3 , A2I2P1
RefSeq Size:	4204
Cytogenetics:	11p11.2
RefSeq ORF:	1779
Synonyms:	HCRY2; PHLL2
Summary:	This gene encodes a flavin adenine dinucleotide-binding protein that is a key component of the circadian core oscillator complex, which regulates the circadian clock. This gene is upregulated by CLOCK/ARNTL heterodimers but then represses this upregulation in a feedback loop using PER/CRY heterodimers to interact with CLOCK/ARNTL. Polymorphisms in this gene have been associated with altered sleep patterns. The encoded protein is widely conserved across plants and animals. Two transcript variants encoding different isoforms have been found for this gene. [provided by RefSeq, Feb 2014]
Protein Families:	Druggable Genome
Protein Pathways:	Circadian rhythm - mammal

Product images:



Coomassie blue staining of purified CRY2 protein (Cat# TP307337). The protein was produced from HEK293T cells transfected with CRY2 cDNA clone (Cat# [RC207337]) using MegaTran 2.0 (Cat# [TT210002]).