

Product datasheet for TP307402

TP53I11 (NM_006034) Human Recombinant Protein


Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human tumor protein p53 inducible protein 11 (TP53I11), 20 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207402 protein sequence Red =Cloning site Green =Tags(s)
	MAAKQPPPLMKKHSQTDLVSRLKTRKILGVGGEDDDGEVHRKISQVLGNEIKFTIREPLGLRVWQFVSA VLFSGIAIMALAFPDQLYDAVFDGAQVTSKTPIRLYGGALLSISLIMWNALYAEKVIIRWTLLTEACYF GVQFLVWTATLAETGLMSLGILLLLVSRLLFVISIYYYYQVGRRPKKA
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	20.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_006025
Locus ID:	9537
UniProt ID:	O14683
RefSeq Size:	3518



[View online »](#)

Cytogenetics: 11p11.2
RefSeq ORF: 567
Synonyms: PIG11
Protein Families: Druggable Genome, Transmembrane

Product images:

Coomassie blue staining of purified TP53I11 protein (Cat# TP307402). The protein was produced from HEK293T cells transfected with TP53I11 cDNA clone (Cat# [RC207402]) using MegaTran 2.0 (Cat# [TT210002]).