

Product datasheet for **TP307432M**

Cytochrome P450 3A5 (CYP3A5) (NM_000777) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human cytochrome P450, family 3, subfamily A, polypeptide 5 (CYP3A5), 100 µg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207432 protein sequence Red =Cloning site Green =Tags(s)
	<p>MDLIPNLAVETWLLLAVSLVLLYLYGTRTHGLFKRLGIPGPTPLPLLGNVLSYRQGLWKFDTECYKKGK MWGTYEGQLPVLAITDPDVIRTVLVKECYSVFTNRRSLGPVGFMKSAISLAEDEEWKRIRLLSPTFTSG KLKEMFPIAQYGDVLRNLRREAEGKPVTLKDIFGAYSMDVITGTSFGVNIDSLNPNQDPFVESTKFF LKFGFLDPLFSLIILFPFLTPVFEALNVSLFPKDTINFLSKSVNRMKKSRLNDKQKHRLDLFLQLMIDSQN SKETESHKALSDLELAAQSIIFIFAGYETTSSVLSFTLYELATHPDVQKQKQKEIDAVLPNKAPPTYDAV VQMEYLDMMVNETLRLFPVAIRLERTCKKDVINGVFIPKGSMMVIPTYALHHDPKYWTEPEEFRPERFS KKKDSIDPYIYTPFGTGPRNCIGMRFALMNMKLALIRVLQNFSPKCKETQIPLKLDLDTQGLLQPEKPIVL KVDSRDGTLTSGE</p> <p>TRTRPLEQKLISEEDLAANDILDYKDDDDKV</p>
Tag:	C-Myc/DDK
Predicted MW:	56.9 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.



[View online »](#)

Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_000768
Locus ID:	1577
UniProt ID:	P20815
RefSeq Size:	1755
Cytogenetics:	7q22.1
RefSeq ORF:	1506
Synonyms:	CP35; CYP11A5; P450PCN3; PCN3
Summary:	This gene encodes a member of the cytochrome P450 superfamily of enzymes. The cytochrome P450 proteins are monooxygenases which catalyze many reactions involved in drug metabolism and synthesis of cholesterol, steroids and other lipids. The encoded protein metabolizes drugs as well as the steroid hormones testosterone and progesterone. This gene is part of a cluster of cytochrome P450 genes on chromosome 7q21.1. Two pseudogenes of this gene have been identified within this cluster on chromosome 7. Expression of this gene is widely variable among populations, and a single nucleotide polymorphism that affects transcript splicing has been associated with susceptibility to hypertension. Alternative splicing results in multiple transcript variants. [provided by RefSeq, Apr 2014]
Protein Families:	Druggable Genome, P450, Transmembrane
Protein Pathways:	Drug metabolism - cytochrome P450, Drug metabolism - other enzymes, Linoleic acid metabolism, Metabolic pathways, Metabolism of xenobiotics by cytochrome P450, Retinol metabolism

Product images:



Coomassie blue staining of purified CYP3A5 protein (Cat# [TP307432]). The protein was produced from HEK293T cells transfected with CYP3A5 cDNA clone (Cat# [RC207432]) using MegaTran 2.0 (Cat# [TT210002]).