

Product datasheet for TP307444


CaMKI gamma (CAMK1G) (NM_020439) Human Recombinant Protein

Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human calcium/calmodulin-dependent protein kinase IG (CAMK1G)
Species:	Human
Expression Host:	HEK293T
Tag:	C-Myc/DDK
Predicted MW:	52.9 kDa
Concentration:	>50 ug/mL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris.HCl, pH 7.3, 100 mM glycine, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_065172
Locus ID:	57172
RefSeq Size:	2507
Cytogenetics:	1q32.2
RefSeq ORF:	1428
Synonyms:	CLICK3; CLICKIII; dj272L16.1; VWS1
Summary:	This gene encodes a protein similar to calcium/calmodulin dependent protein kinase, however, its exact function is not known. [provided by RefSeq, Jul 2008]
Protein Families:	Druggable Genome, Protein Kinase



[View online »](#)

Product images:

Coomassie blue staining of purified CAMK1G protein (Cat# TP307444). The protein was produced from HEK293T cells transfected with CAMK1G cDNA clone (Cat# [RC207444]) using MegaTran 2.0 (Cat# [TT210002]).