

Product datasheet for **TP307463L**

DEXI (NM_014015) Human Recombinant Protein

Product data:


Product Type:	Recombinant Proteins
Description:	Recombinant protein of human dexamethasone-induced transcript (DEXI), 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207463 protein sequence Red =Cloning site Green =Tags(s)
	MLGARVA AHLDALG PLVYPVPPLLPSMFYVGLFFVNV LILYYAFLMEYIVLNVGLVFLPEDMDQALVDL GVLSDPGSGLYDADSELDFDAYLE
	TRTRPLEQKLISEEDLAANDILDYKDDDDKV
Tag:	C-Myc/DDK
Predicted MW:	10.2 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	NP_054734
Locus ID:	28955
UniProt ID:	O95424
RefSeq Size:	1580
Cytogenetics:	16p13.13



[View online »](#)

RefSeq ORF: 285
Synonyms: MYLE
Protein Families: Transmembrane

Product images:



Coomassie blue staining of purified DEXI protein (Cat# [TP307463]). The protein was produced from HEK293T cells transfected with DEXI cDNA clone (Cat# [RC207463]) using MegaTran 2.0 (Cat# [TT210002]).