

## Product datasheet for TP307464L

### KCNIP4 (NM\_001035004) Human Recombinant Protein

#### Product data:

Product Type:	Recombinant Proteins
Description:	Recombinant protein of human Kv channel interacting protein 4 (KCNIP4), transcript variant 6, 1 mg
Species:	Human
Expression Host:	HEK293T
Expression cDNA Clone or AA Sequence:	>RC207464 protein sequence Red=Cloning site Green=Tags(s)

MNVRRVESISAQLEEASSTGGFLYAQNSTKRSIKERLMKLLPCSAAKTSSPAIQNSVEDELEMATVRHRP  
EALELLEAQSKFTKKELQILYRGFKNECPSGVNEETFKEIYSQFFPQGDSTTYAHFLFNAFDTDHNGAV  
SFEDFIKGLSILLRGTVQEKLNWAFNLYDINKDGYITKEEMLDIMKAIYDMMGKCTYPVLKEDAPRQHVE  
TFFQKMDKNKDGWVTIDEFIESCQKDENIMRSMQLFENVI

TRTRPLEQKLISEEDLAANDILDYKDDDDKV

Tag:	C-Myc/DDK
Predicted MW:	21.8 kDa
Concentration:	>0.05 µg/µL as determined by microplate BCA method
Purity:	> 80% as determined by SDS-PAGE and Coomassie blue staining
Buffer:	25 mM Tris-HCl, 100 mM glycine, pH 7.3, 10% glycerol
Preparation:	Recombinant protein was captured through anti-DDK affinity column followed by conventional chromatography steps.
Note:	For testing in cell culture applications, please filter before use. Note that you may experience some loss of protein during the filtration process.
Storage:	Store at -80°C.
Stability:	Stable for 12 months from the date of receipt of the product under proper storage and handling conditions. Avoid repeated freeze-thaw cycles.
RefSeq:	<a href="#">NP_001030176</a>
Locus ID:	80333



[View online »](#)

UniProt ID: [Q6PIL6](#), [A8K4N3](#)

RefSeq Size: 2345

Cytogenetics: 4p15.31-p15.2


RefSeq ORF: 753

Synonyms: CALP; KCHIP4

**Summary:** This gene encodes a member of the family of voltage-gated potassium (Kv) channel-interacting proteins (KCNIPs), which belong to the recoverin branch of the EF-hand superfamily. Members of the KCNIP family are small calcium binding proteins. They all have EF-hand-like domains, and differ from each other in the N-terminus. They are integral subunit components of native Kv4 channel complexes. They may regulate A-type currents, and hence neuronal excitability, in response to changes in intracellular calcium. This protein member also interacts with presenilin. Multiple alternatively spliced transcript variants encoding distinct isoforms have been identified for this gene. [provided by RefSeq, Jul 2008]

**Protein Families:** Druggable Genome, Ion Channels: Other

### Product images:



Coomassie blue staining of purified KCNIP4 protein (Cat# [TP307464]). The protein was produced from HEK293T cells transfected with KCNIP4 cDNA clone (Cat# [RC207464]) using MegaTran 2.0 (Cat# [TT210002]).